

# Intuitive imaging, best results.

Camera Module CM 900®



Digital Slit  
Lamp Imaging

 **HAAG-STREIT**  
INTERNATIONAL  
Precision by Tradition

# The intelligent choice for optimal imaging.

- Easy intuitive imaging
- Still imaging and high resolution video<sup>1)</sup>
- Freeze Technology – high speed capture and data transfer
- Easy Touch – precise image control
- Wide dynamic range, high sensitivity
- History Trigger – never miss the perfect image
- Compact

<sup>1)</sup> Video capture is only available with EyeCap® SL

## From the leader in slit lamp imaging



Slitlamp BD 900®

EASY TOUCH (optional)



**Control the image:** Easy touch more than just a trigger, allowing simple, ergonomic management of the camera settings.

Slitlamp BP 900®  
(In combination with Video Adapter)

Camera Module CM 900®, exceptional imaging solution for the BD 900®, BP 900® and BMV 900®<sup>2)</sup> slit lamp. Combination of the CM 900® with EyeCap™ or EyeSuite™ Imaging completes the perfect imaging system.


<sup>2)</sup> BMV discontinued

## FREEZE TECHNOLOGY



**Freeze technology:** Captures at the precise moment you press the trigger.

## HISTORY TRIGGER\*



**Rewind time and freeze the moment:** The history trigger not only captures one image but records the last few seconds in real time and lets you return to select the optimum image.

\* Patent pending

**HAAG-STREIT AG**  
Gartenstadtstrasse 10  
CH-3098 Koeniz  
Tel. +41 31 978 01 11  
Fax +41 31 978 02 82  
info@haag-streit.ch  
www.haag-streit.com

 **HAAG-STREIT  
INTERNATIONAL**  
Precision by Tradition