

# Customization

PST's instrument trays can be customized with our unique silicone bar positioning system

PST instrument trays can be customized to adapt to a wide variety of individualized configurations. Our uniquely engineered GridLok System™ allows for easy positioning of the special Plug'nSert™ silicone insert bars. PST silicone products are molded from medical grade silicone. PST trays are manufactured with state-of-the-art durable medical grade plastics. PST trays are engineered with a "pebble-like" interior surface to insure proper penetration of steam and other sterilizing agents under the mat.



The PST GridLok System provides maximum flexibility for custom applications. Plug'nSert silicone bars engineered for use with the GridLok System are designed to be inserted onto the tray surface. They simply plug in with a minimum of effort.



For more information or to order  
CALL: 1-800-PST-TRAY

**Customization:**

PST trays can be customized for individual needs.

**Protection:**

PST Trays protect your instruments from costly damage.

**Quality/Durability:**

PST Trays are manufactured with strict quality standards and are engineered to last.

**Versatility:**


PST Trays can be used with a wide variety of instruments and can be individually customized.

**Cost Savings:**


PST's unique Gridlock System™ allows for the use of one tray for many applications.

**Sterility Assurance:**

The raised mat assures proper penetration of sterilizing agents between the mat and the tray surface.



PST Trays are engineered with a unique “pebble-like” interior surface. This raises the mat to assure proper penetration of the sterilizing agent between the mat and the tray surface. Independent laboratory tests have shown that the space between the mat and the surface is effectively penetrated.



Our latch is designed to last! Compare the quality and manufacturing standards of PST engineered trays and see the difference. Mold flow analysis gives PST products the strength and durability to satisfy your requirements.

Surgical Instrument Trays By PST			
Part #	Description	Dimensions W x L x D	
		Inches	Centimeters
<b>MINI TRAYS</b>			
1300A	Mini Tray - Base, Lid, & Mat	1.5x2.7x0.75	3.8x6.8x1.9
1310A	Mini Tray - Base (deep), Lid, & Mat	1.5x2.7x1.0	3.8x6.8x2.5
1320A	Mini Tray - Base (deep), Lid, & Mat	1.5x2.7x1.25	3.8x6.8x3.2
<b>LARGE MICROSURGICAL INSTRUMENT TRAYS</b>			
1500A	Instrument Tray - Base, Lid, & Mat	10x15x0.75	25.4x38.1x1.9
1520A	Instrument Tray - Base, Insert Tray, Lid, & 2 Mats	10x15x1.5	25.4x38.1x3.8
1522A	Instrument Tray - Base (deep), Lid, & Mat	10x15x1.5	25.4x38.1x3.8
<b>PHACO TRAYS</b>			
2341A	Phaco Tray - Base, Lid, & 2 Bars / 1 Slot	4.0x7.5x1.5	10.1x19x3.8
2343A	Phaco Tray - Base, Lid, & 2 Bars / 3 Slot	4.0x7.5x1.5	10.1x19x3.8
2361A	Phaco Tray - Base, Lid, & 2 Bars / 1 Slot	6.0x10x1.5	15.2x25.4x3.8
2363A	Phaco Tray - Base, Lid, & 2 Bars / 3 Slot	6.0x10x1.5	15.2x25.4x3.8
<b>DIAMOND KNIFE TRAYS</b>			
2501A	Knife Tray, 1 Slot - Base, Lid, & 2 Bars	2.5x6.0x0.75	6.4x15.2x1.9
2502A	Knife Tray, 2 Slot - Base, Lid, & 2 Bars	2.5x6.0x0.75	6.4x15.2x1.9
2503A	Knife Tray, 3 Slot - Base, Lid, & 2 Bars	2.5x6.0x0.75	6.4x15.2x1.9
2505A	Knife Tray, 5 Slot - Base, Lid, & 2 Bars	4.0x6.5x0.75	10.2x16.5x1.9
2571A	Knife Tray, 1 Slot - Base, Lid, & 2 Bars	2.5x7.5x0.75	6.4x19x1.9
2572A	Knife Tray, 2 Slot - Base, Lid, & 2 Bars	2.5x7.5x0.75	6.4x19x1.9
2573A	Knife Tray, 3 Slot - Base, Lid, & 2 Bars	2.5x7.5x0.75	6.4x19x1.9
<b>MICROSURGICAL INSTRUMENT TRAYS</b>			
2600A	Micro Tray - Base, Lid, & Mat	2.5x6.0x0.75	6.4x15.2x1.9
2610A	Micro Tray - Base, Lid, & Mat	2.5x6.0x1.25	6.4x15.2x3.2
2700A	Micro Tray - Base, Lid, & Mat	2.5x7.5x0.75	6.4x19x1.9
4600A	Micro Tray - Base, Lid, & Mat	4.0x6.5x0.75	10.2x16.2x1.9
4700A	Micro Tray - Base, Lid, & Mat	4.0x7.5x0.75	10.1x19x1.9
4710A	Micro Tray - Base, Lid, & Mat	4.0x7.5x1.5	10.1x19x3.8
6100A	Micro Tray - Base, Lid, & Mat	6.0x10x0.75	15.2x25.4x1.9
6120A	Micro Tray - Base, Insert Tray, Lid, & 2 Mats	6.0x10x1.5	15.2x25.4x3.8
6124A	Micro Tray - Base (deep), Lid, & Mat	6.0x10x1.5	15.2x25.4x3.8

PST Corp. warrants its trays and components for one year from the date of purchase against failure due to manufacturing defects. PST will replace defective products at no charge to the customer upon prepaid return of the product to PST.

2600X 1C2.5M CE



For more information or to order  
CALL: 1-800-PST-TRAY